

DATASHEET

CD7 Rabbit Monoclonal Antibody(C682)

CAT. NO. AMA00294

KEY FEATURES

| | | | |
|---------------|---|---------------|--------------|
| Target | CD7 | Source / Host | Rabbit |
| Reactivity | Human | Clonality | Monoclonal |
| Applications | IHC | Conjugation | Unconjugated |
| Form / Buffer | Liquid in PBS, pH 7.4, containing 50% glycerol, 0.05% BSA and 0.01% sodium azide. | Storage | at-20°C |

BACKGROUND

Transmembrane glycoprotein expressed by T-cells and natural killer (NK) cells and their precursors cells and their precursors . Plays a costimulatory role in T-cell activation upon binding to its ligand K12/SECTM1 . In turn, mediates the production of cytokines such as IL-2 . On resting NK-cells, CD7 activation results in a significant induction of interferon-gamma levels .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

| | |
|-----|---------------|
| IHC | 1:100 - 1:500 |
|-----|---------------|

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


| | |
|-------------------|--|
| Description | Recombinant rabbit monoclonal antibody to CD7 |
| Specificity | Recognizes endogenous levels of CD7 protein |
| Antibody Type | Primary antibody, Recombinant |
| Immunogen | Recombinant fusion protein of human CD7. The exact sequence is proprietary. |
| Purification | The antibody was purified by immunogen affinity chromatography. |
| Form/Buffer | Liquid in PBS, pH 7.4, containing 50% glycerol, 0.05% BSA and 0.01% sodium azide. |
| Alternative Names | T-cell antigen CD7; GP40; T-cell leukemia antigen; T-cell surface antigen Leu-9; TP41; CD7 |
| Gene Symbol | CD7 |
| Entrez Gene | 924(Human) |
| SwissProt | P09564(Human) |

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**CD7 Rabbit Monoclonal Antibody(C682)**

CAT. NO. AMA00294

DATA

Immunohistochemical analysis of CD7 staining in human tonsil formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.