

## DATASHEET

# PIP Rabbit Monoclonal Antibody(C706)

CAT. NO. AMA00318

### KEY FEATURES

Target	PIP	Source / Host	Rabbit
Reactivity	Human	Clonality	Monoclonal
Applications	IHC	Conjugation	Unconjugated
Form / Buffer	Liquid in PBS, pH 7.4, containing 50% glycerol, 0.05% BSA and 0.01% sodium azide.	Storage	at-20°C

### BACKGROUND

Prolactin-inducible protein (PIP), also known as Gross cystic disease fluid protein 15, GCDFP-15, Prolactin-induced protein. Belongs to the PIP family Subunit: Monomer. Interacts with AZGP1 Subcellular location: Secreted Tissue specificity: Expressed in pathological conditions of the mammary gland and in several exocrine tissues, such as the lacrimal, salivary, and sweat glands Induction: By prolactin and androgen; inhibited by estrogen

### APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

IHC	1:50 - 1:300
-----	--------------

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

### PRODUCT OVERVIEW


Description	Recombinant rabbit monoclonal antibody to PIP
Specificity	Recognizes endogenous levels of PIP protein
Antibody Type	Primary antibody, Recombinant
Immunogen	Recombinant fusion protein of human PIP. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Form/Buffer	Liquid in PBS, pH 7.4, containing 50% glycerol, 0.05% BSA and 0.01% sodium azide.
Alternative Names	GCDFP15; GPIP4; Prolactin-inducible protein; Gross cystic disease fluid protein 15; GCDFP-15; Prolactin-induced protein; Secretory actin-binding protein; SABP; gp17
Gene Symbol	PIP
Entrez Gene	5304(Human)
SwissProt	P12273(Human)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****PIP Rabbit Monoclonal Antibody(C706)**

CAT. NO. AMA00318

**DATA**

Immunohistochemical analysis of PIP staining in human breast cancer formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.