

DATASHEET

MMP14 Rabbit Monoclonal Antibody(C1548)

CAT. NO. AMA01160

KEY FEATURES

Target	MMP14	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Monoclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in PBS, pH 7.4, containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.	Storage	at-20°C

BACKGROUND

Endopeptidase that degrades various components of the extracellular matrix such as collagen . Essential for pericellular collagenolysis and modeling of skeletal and extraskelatal connective tissues during development . Activates progelatinase A/MMP2, thereby acting as a positive regulator of cell growth and migration . Involved in the formation of the fibrovascular tissues in association with pro-MMP2 . May be involved in actin cytoskeleton reorganization by cleaving PTK7 . Acts as a regulator of Notch signaling by mediating cleavage and inhibition of DLL1 . Cleaves ADGRB1 to release vasculostatin-40 which inhibits angiogenesis . Acts as a negative regulator of the GDF15-GFRAL aversive response by mediating cleavage and inactivation of GFRAL .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Recombinant rabbit monoclonal antibody to MMP14
Specificity	Recognizes endogenous levels of MMP14 protein
Antibody Type	Primary antibody, Recombinant
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human MMP14. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 65 kD; Observed: 65 kD
Form/Buffer	Liquid in PBS, pH 7.4, containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.
Alternative Names	Matrix metalloproteinase-14; MMP-14; MMP-X1; Membrane-type matrix metalloproteinase 1; MT-MMP 1; MTMMP1; Membrane-type-1 matrix metalloproteinase; MT1-MMP; MT1MMP
Gene Symbol	MMP14
Entrez Gene	4323(Human); 17387(Mouse); 81707(Rat)
SwissProt	P50281(Human); P53690(Mouse); Q10739(Rat)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**MMP14 Rabbit Monoclonal Antibody(C1548)**

CAT. NO. AMA01160

DATA

Western blot analysis of MMP14 expression in EC9706 (A), HeLa (B), THP1 (C), mouse lung (D), mouse kidney (E), rat brain (F) whole cell lysates. (Predicted band size: 65 kD; Observed band size: 65 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.