

**DATASHEET**

**TRMT2A Rabbit Monoclonal Antibody(C1660)**

CAT. NO. AMA01272

**KEY FEATURES**

Target	TRMT2A	Source / Host	Rabbit
Reactivity	Human	Clonality	Monoclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in PBS, pH 7.4, containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.	Storage	at-20°C

**BACKGROUND**

S-adenosyl-L-methionine-dependent methyltransferase that catalyzes the formation of 5-methyl-uridine in tRNAs and some mRNAs . Mainly catalyzes the methylation of uridine at position 54 (m5U54) in cytosolic tRNAs . Also able to mediate the formation of 5-methyl-uridine in some mRNAs .

**APPLICATION**

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
----	----------------

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

**PRODUCT OVERVIEW**


Description	Recombinant rabbit monoclonal antibody to TRMT2A
Specificity	Recognizes endogenous levels of TRMT2A protein
Antibody Type	Primary antibody, Recombinant
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human TRMT2A protein. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 68 kD; Observed: 75 kD
Form/Buffer	Liquid in PBS, pH 7.4, containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.
Alternative Names	HTF9C; tRNA (uracil-5-)-methyltransferase homolog A; HpaII tiny fragments locus 9c protein
Gene Symbol	TRMT2A
Entrez Gene	27037(Human)
SwissProt	Q8IZ69(Human)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****TRMT2A Rabbit Monoclonal Antibody(C1660)**

CAT. NO. AMA01272

**DATA**

Western blot analysis of TRMT2A expression in HEK293T (A), MDAMB231 (B) whole cell lysates. (Predicted band size: 68 kD; Observed band size: 75 kD)

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.