

DATASHEET

COASY Rabbit Monoclonal Antibody(C1673)

CAT. NO. AMA01285

KEY FEATURES

Target	COASY	Source / Host	Rabbit
Reactivity	Human, Rat	Clonality	Monoclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in PBS, pH 7.4, containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.	Storage	at-20°C

BACKGROUND

Bifunctional enzyme that catalyzes the fourth step of the coenzyme A biosynthetic pathway, the adenylation of 4'-phosphopantetheine, and the fifth step, the phosphorylation of dephospho-CoA to CoA.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Recombinant rabbit monoclonal antibody to COASY
Specificity	Recognizes endogenous levels of COASY protein
Antibody Type	Primary antibody, Recombinant
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human COASY protein. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 62 kD; Observed: 62 kD
Form/Buffer	Liquid in PBS, pH 7.4, containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.
Alternative Names	Bifunctional coenzyme A synthase; CoA synthase; NBP; POV-2
Gene Symbol	COASY
Entrez Gene	80347(Human)
SwissProt	Q13057(Human)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**COASY Rabbit Monoclonal Antibody(C1673)**

CAT. NO. AMA01285

DATA

Western blot analysis of COASY expression in rat kidney (A), HEK293T (B), MDAMB231 (C) whole cell lysates. (Predicted band size: 62 kD; Observed band size: 62 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.