

DATASHEET

DNAJC12 Rabbit Monoclonal Antibody(C1697)

CAT. NO. AMA01309

KEY FEATURES

Target	DNAJC12	Source / Host	Rabbit
Reactivity	Human, Mouse	Clonality	Monoclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in PBS, pH 7.4, containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.	Storage	at-20°C

BACKGROUND

Probable co-chaperone that participates in the proper folding of bipterin-dependent aromatic amino acid hydroxylases, which include phenylalanine-4-hydroxylase (PAH), tyrosine 3-monooxygenase (TH) and peripheral and neuronal tryptophan hydroxylases (TPH1 and TPH2).

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Recombinant rabbit monoclonal antibody to DNAJC12
Specificity	Recognizes endogenous levels of DNAJC12 protein
Antibody Type	Primary antibody, Recombinant
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human DNAJC12 protein. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 23 kD; Observed: 23 kD
Form/Buffer	Liquid in PBS, pH 7.4, containing 50% glycerol, 0.2% BSA and 0.01% sodium azide.
Alternative Names	JDP1; Dnaj homolog subfamily C member 12; J domain-containing protein 1
Gene Symbol	DNAJC12
Entrez Gene	56521(Human); 30045(Mouse)
SwissProt	Q9UKB3(Human); Q9R022(Mouse)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**DNAJC12 Rabbit Monoclonal Antibody(C1697)**

CAT. NO. AMA01309

DATA

Western blot analysis of DNAJC12 expression in mouse kidney (A) whole cell lysates.
(Predicted band size: 23 kD; Observed band size: 23 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.