

DATASHEET

TMEM256 Rabbit Monoclonal Antibody(C1985)

CAT. NO. AMA01597

KEY FEATURES

Target	TMEM256	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Monoclonal
Applications	WB, IP, FC	Conjugation	Unconjugated
Form / Buffer	Liquid in PBS, pH 7.3, 50% glycerol, 0.05% BSA, and 0.05% Proclin300.	Storage	at-20°C

BACKGROUND

Transmembrane protein 256 (TMEM256). Belongs to the TMEM256 family Subcellular location: Membrane

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
IP	1:10 - 1:50
FC	1:50 - 1:100

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Recombinant rabbit monoclonal antibody to TMEM256
Specificity	Recognizes endogenous levels of TMEM256 protein
Antibody Type	Primary antibody, Recombinant
Immunogen	Recombinant fusion protein of human TMEM256. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 11 kD; Observed: 12 kD
Form/Buffer	Liquid in PBS, pH 7.3, 50% glycerol, 0.05% BSA, and 0.05% Proclin300.
Alternative Names	C17orf61; Transmembrane protein 256
Gene Symbol	TMEM256
Entrez Gene	254863(Human); 69186(Mouse)
SwissProt	Q8N2U0(Human); Q5F285(Mouse)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**TMEM256 Rabbit Monoclonal Antibody(C1985)**

CAT. NO. AMA01597

DATA

Western blot analysis of TMEM256 expression in HeLa (A) whole cell lysates. (Predicted band size: 11 kD; Observed band size: 12 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.