

DATASHEET

GFP-tag Mouse Monoclonal Antibody(C2041)

CAT. NO. AMA01653

KEY FEATURES

Target	GFP-tag	Source / Host	Mouse
Reactivity	PlantSheep	Clonality	Monoclonal
Applications	WB, IP	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at-20°C

BACKGROUND

Energy-transfer acceptor. Its role is to transduce the blue chemiluminescence of the protein aequorin into green fluorescent light by energy transfer. Fluoresces in vivo upon receiving energy from the Ca(2+)-activated photoprotein aequorin.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:2000 - 1:5000
IP	1:100 - 1:200

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Mouse monoclonal antibody to GFP-tag
Specificity	Recognizes GFP and GFP tag fusion proteins.
Antibody Type	Primary antibody,Tag
Immunogen	Recombinant protein corresponding to GFP-tag.
Purification	Affinity chromatography
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
SwissProt	P42212(Human)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.


*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**GFP-tag Mouse Monoclonal Antibody(C2041)**

CAT. NO. AMA01653

DATA

Western blot analysis of over-expressed GFP-tagged protein in 293T cell lysate.



Immunoprecipitation of GFP-tagged protein from HEK293T cells transfected with vector overexpressing GFP tag, using Anti-GFP-tag Antibody.

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.