

**DATASHEET**

**Biotin-conjugated Hepatitis B Virus E Antigen Mouse Monoclonal Antibody(C2467)**

CAT. NO. AMA02079

**KEY FEATURES**

Target	Hepatitis B Virus E Antigen	Source / Host	Mouse
Reactivity	Hepatitis B virus	Clonality	Monoclonal
Applications	ELISA	Conjugation	Biotin
Form / Buffer	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.	Storage	at-20°C

**BACKGROUND**

Hepatitis B e antigen (HBeAg) is a secreted form of the HBV core protein, considered a marker of active viral replication and high infectivity. HBeAg seroconversion (HBeAg loss with anti-HBe production) indicates reduced viral activity and a more favorable prognosis. Anti-HBeAg antibodies are used in serological diagnosis and monitoring of chronic HBV infection.

**APPLICATION**

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

ELISA	Use at an assay dependent dilution.
-------	-------------------------------------

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

**PRODUCT OVERVIEW**

Description	Mouse monoclonal antibody to Hepatitis B Virus E Antigen-Biotin labeled
Specificity	Reacts only with Hepatitis B Virus E Antigen.
Antibody Type	Primary antibody
Immunogen	Recombinant Hepatitis B Virus E Antigen
Purification	Immunogen affinity purified
Form/Buffer	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Alternative Names	Core antigen; Core protein; HBe antigen; HBeAg; HBVgp4; Pre C C; Precore protein; Precore/core; Precore/core ORF; Precore/core protei

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****Biotin-conjugated Hepatitis B Virus E Antigen Mouse Monoclonal Antibody(C2467)**

CAT. NO. AMA02079

**| DATA****| STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**| NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.