

**DATASHEET**

**Histone Deacetylase 4 Mouse Monoclonal Antibody(C2735)**

CAT. NO. AMA02347

**KEY FEATURES**

Target	Histone Deacetylase 4	Source / Host	Mouse
Reactivity	Human	Clonality	Monoclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Mouse IgG1. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.	Storage	at-20°C

**BACKGROUND**

Responsible for the deacetylation of lysine residues on the N-terminal part of the core histones (H2A, H2B, H3 and H4). Histone deacetylation gives a tag for epigenetic repression and plays an important role in transcriptional regulation, cell cycle progression and developmental events. Histone deacetylases act via the formation of large multiprotein complexes. Involved in muscle maturation via its interaction with the myocyte enhancer factors such as MEF2A, MEF2C and MEF2D. Involved in the MTA1-mediated epigenetic regulation of ESR1 expression in breast cancer. Deacetylates HSPA1A and HSPA1B at 'Lys-77' leading to their preferential binding to co-chaperone STUB1 .

**APPLICATION**

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
----	----------------

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

**PRODUCT OVERVIEW**


Description	Mouse monoclonal to Histone Deacetylase 4
Specificity	Recognizes endogenous levels of Histone Deacetylase 4 protein
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human Histone Deacetylase 4 expressed in E. Coli
Purification	This antibody is purified through a protein G column.
Molecular Weight	Predicted: 119 kD; Observed: 110 kD kD
Form/Buffer	Mouse IgG1. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	KIAA0288; Histone deacetylase 4; HD4
Gene Symbol	HDAC4
Entrez Gene	9759(Human)
SwissProt	P56524(Human)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****Histone Deacetylase 4 Mouse Monoclonal Antibody(C2735)**

CAT. NO. AMA02347

**DATA**

Western blot analysis of Histone Deacetylase 4 expression in HeLa (A), Jurkat (B) whole cell lysates. (Predicted band size: 119 kD; Observed band size: 110 kD kD)

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.