

DATASHEET

CD158e Mouse Monoclonal Antibody(C2796)

CAT. NO. AMA02408

KEY FEATURES

Target	CD158e	Source / Host	Mouse
Reactivity	Human, Rat	Clonality	Monoclonal
Applications	WB, FC	Conjugation	Unconjugated
Form / Buffer	Mouse IgG2b. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.	Storage	at-20°C

BACKGROUND

Receptor on natural killer (NK) cells for HLA Bw4 allele. Inhibits the activity of NK cells thus preventing cell lysis.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
FC	1:100 - 1:200

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

Description	Mouse monoclonal to CD158e
Specificity	Recognizes endogenous levels of CD158e protein
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human CD158e expressed in HEK293
Purification	This antibody is purified through a protein G column.
Molecular Weight	Predicted: 49 kD; Observed: 49 kD
Form/Buffer	Mouse IgG2b. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	CD158E; NKAT3; NKB1; Killer cell immunoglobulin-like receptor 3DL1; CD158 antigen-like family member E; HLA-BW4-specific, inhibitory NK cell receptor; MHC class I NK cell receptor; Natural killer-associated transcript 3; NKAT-3; p70 natural killer cell receptor clones CL-2/CL-11; p70 NK receptor CL-2/CL-11; CD158e
Gene Symbol	KIR3DL1
Entrez Gene	3811(Human); 353253(Rat)
SwissProt	P43629(Human); P83556(Rat)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.


*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET


CD158e Mouse Monoclonal Antibody(C2796)

CAT. NO. AMA02408

DATA



Western blot analysis of CD158e expression in A431 (A), Raji (B), SPCA1 (C), K562 (D), HEK293 (E), U937 (F), C6 (G) whole cell lysates. (Predicted band size: 49 kD; Observed band size: 49 kD)



Flow cytometric analysis of Jurkat cells using Anti-CD158e Antibody (green) and negative control (red).

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.