

DATASHEET

RAD50 Mouse Monoclonal Antibody(C2931)

CAT. NO. AMA02543

KEY FEATURES

Target	RAD50	Source / Host	Mouse
Reactivity	Human, Rat	Clonality	Monoclonal
Applications	WB, IHC, FC	Conjugation	Unconjugated
Form / Buffer	Mouse IgG2a. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.	Storage	at-20°C

BACKGROUND

Component of the MRN complex, which plays a central role in double-strand break (DSB) repair, DNA recombination, maintenance of telomere integrity and meiosis repair, DNA recombination, maintenance of telomere integrity and meiosis . The MRN complex is involved in the repair of DNA double-strand breaks (DSBs) via homologous recombination (HR), an error-free mechanism which primarily occurs during S and G2 phases . The complex (1) mediates the end resection of damaged DNA, which generates proper single-stranded DNA, a key initial steps in HR, and is (2) required for the recruitment of other repair factors and efficient activation of ATM and ATR upon DNA damage .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
IHC	1:100 - 1:500
FC	1:100 - 1:200

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

Description	Mouse monoclonal to RAD50
Specificity	Recognizes endogenous levels of RAD50 protein
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human RAD50 expressed in E. Coli
Purification	This antibody is purified through a protein G column.
Molecular Weight	Predicted: 154 kD; Observed: 135 kD kD
Form/Buffer	Mouse IgG2a. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	DNA repair protein RAD50; hRAD50
Gene Symbol	RAD50
Entrez Gene	10111(Human); 64012(Rat)
SwissProt	Q92878(Human); Q9JIL8(Rat)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.


*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET


RAD50 Mouse Monoclonal Antibody(C2931)

CAT. NO. AMA02543


DATA



Western blot analysis of RAD50 expression in C6 (A), HepG2 (B) whole cell lysates. (Predicted band size: 154 kD; Observed band size: 135 kD)



Immunohistochemical analysis of RAD50 staining in human esophageal cancer formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.



Flow cytometric analysis of HeLa cells using Anti-RAD50 Antibody (green) and negative control (red).

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.