

DATASHEET

FKBP4 Rabbit Monoclonal Antibody(C3149)

CAT. NO. AMA02761

KEY FEATURES

Target	FKBP4	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat, Hamster	Clonality	Monoclonal
Applications	WB, IF/ICC	Conjugation	Unconjugated
Form / Buffer	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% BSA.	Storage	at-20°C

BACKGROUND

Immunophilin protein with PPIase and co-chaperone activities. Component of steroid receptors heterocomplexes through interaction with heat-shock protein 90 (HSP90). May play a role in the intracellular trafficking of heterooligomeric forms of steroid hormone receptors between cytoplasm and nuclear compartments. The isomerase activity controls neuronal growth cones via regulation of TRPC1 channel opening. Also acts as a regulator of microtubule dynamics by inhibiting MAPT/TAU ability to promote microtubule assembly. May have a protective role against oxidative stress in mitochondria.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
IF/ICC	1:50 - 1:100

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

Description	Rabbit monoclonal antibody to FKBP4
Specificity	Recognizes endogenous levels of FKBP4 protein.
Antibody Type	Primary antibody
Immunogen	A synthetic peptide of human FKBP52
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 52 kD; Observed: 52 kD
Form/Buffer	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% BSA.
Alternative Names	FKBP52; Peptidyl-prolyl cis-trans isomerase FKBP4; PPIase FKBP4; 51 kDa FK506-binding protein; FKBP51; 52 kDa FK506-binding protein; 52 kDa FKBP; FKBP-52; 59 kDa immunophilin; p59; FK506-binding protein 4; FKBP-4; FKBP59; HSP-binding immunophilin; HBI; Immunophilin FKBP52; Rotamase
Gene Symbol	FKBP4
Entrez Gene	2288(Human); 14228(Mouse); 260321(Rat)
SwissProt	Q02790(Human); P30416(Mouse); Q9QVC8(Rat)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.


*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET


FKBP4 Rabbit Monoclonal Antibody(C3149)

CAT. NO. AMA02761

DATA



Western blot analysis of FKBP4 expression in Raw264.7 (A), Jurkat (B), C6 (C), CHOK1 (D), HeLa (E) whole cell lysates. (Predicted band size: 52 kD; Observed band size: 52 kD)



Immunofluorescent analysis of FKBP4 staining in HeLa cells. Formalin-fixed cells were permeabilized with 0.1% Triton X-100 in TBS for 5-10 minutes and blocked with 3% BSA-PBS for 30 minutes at room temperature. Cells were probed with the primary antibody in 3% BSA-PBS and incubated overnight at 4 °C in a humidified chamber. Cells were washed with PBST and incubated with a AREX® Fluor 488 -conjugated secondary antibody (green) in PBS at room temperature in the dark.

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.