

DATASHEET

GSDMD-CT Rabbit Monoclonal Antibody(C3159)

CAT. NO. AMA02771

KEY FEATURES

Target	GSDMD-CT	Source / Host	Rabbit
Reactivity	Human	Clonality	Monoclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% BSA.	Storage	at-20°C

BACKGROUND

Precursor of a pore-forming protein that plays a key role in host defense against pathogen infection and danger signals . This form constitutes the precursor of the pore-forming protein: upon cleavage, the released N-terminal moiety (Gasdermin-D, N-terminal) binds to membranes and forms pores, triggering pyroptosis .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Rabbit monoclonal antibody to GSDMD-CT
Specificity	Recognizes endogenous levels of GSDMD-CT protein.
Antibody Type	Primary antibody
Immunogen	Recombinant protein of human GSDMD
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 53 kD; Observed: 53 kD
Form/Buffer	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% BSA.
Alternative Names	DFNA5L; GSDMDC1; Gasdermin-D; Gasdermin domain-containing protein 1
Gene Symbol	GSDMD
Entrez Gene	79792(Human)
SwissProt	P57764(Human)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**GSDMD-CT Rabbit Monoclonal Antibody(C3159)**

CAT. NO. AMA02771

DATA

Western blot analysis of GSDMD-CT expression in Jurkat (A) whole cell lysates. (Predicted band size: 53 kD; Observed band size: 53 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.