

DATASHEET

TRAF6 Rabbit Monoclonal Antibody(C3202)

CAT. NO. AMA02814

KEY FEATURES

Target	TRAF6	Source / Host	Rabbit
Reactivity	Human, Rat	Clonality	Monoclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% BSA.	Storage	at-20°C

BACKGROUND

E3 ubiquitin ligase that, together with UBE2N and UBE2V1, mediates the synthesis of 'Lys-63'-linked-polyubiquitin chains conjugated to proteins, such as ECSIT, IKBKG, IRAK1, AKT1 and AKT2 . Also mediates ubiquitination of free/unanchored polyubiquitin chain that leads to MAP3K7 activation . Leads to the activation of NF-kappa-B and JUN . Seems to also play a role in dendritic cells (DCs) maturation and/or activation . Represses c-Myb-mediated transactivation, in B-lymphocytes . Adapter protein that seems to play a role in signal transduction initiated via TNF receptor, IL-1 receptor and IL-17 receptor . Regulates osteoclast differentiation by mediating the activation of adapter protein complex 1 (AP-1) and NF-kappa-B, in response to RANK-L stimulation .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Rabbit monoclonal antibody to TRAF6
Specificity	Recognizes endogenous levels of TRAF6 protein.
Antibody Type	Primary antibody
Immunogen	A synthetic peptide of human TRAF6
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 60 kD; Observed: 60 kD
Form/Buffer	Liquid in 50mM Tris-Glycine (pH 7.4), 0.15M NaCl, 50% Glycerol, 0.01% Sodium azide and 0.05% BSA.
Alternative Names	RNF85; TNF receptor-associated factor 6; E3 ubiquitin-protein ligase TRAF6; Interleukin-1 signal transducer; RING finger protein 85
Gene Symbol	TRAF6
Entrez Gene	7189(Human); 311245(Rat)
SwissProt	Q9Y4K3(Human); B5DF45(Rat)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**TRAF6 Rabbit Monoclonal Antibody(C3202)**

CAT. NO. AMA02814

DATA

Western blot analysis of TRAF6 expression in K562 (A), rat brain (B), C6 (C), HeLa (D) whole cell lysates. (Predicted band size: 60 kD; Observed band size: 60 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.