

DATASHEET

hnRNP U Mouse Monoclonal Antibody(C3326)

CAT. NO. AMA02938

KEY FEATURES

Target	hnRNP U	Source / Host	Mouse
Reactivity	Human	Clonality	Monoclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.	Storage	at-20°C

BACKGROUND

DNA- and RNA-binding protein involved in several cellular processes such as nuclear chromatin organization, telomere-length regulation, transcription, mRNA alternative splicing and stability, Xist-mediated transcriptional silencing and mitotic cell progression . Plays a role in the regulation of interphase large-scale gene-rich chromatin organization through chromatin-associated RNAs (caRNAs) in a transcription-dependent manner, and thereby maintains genomic stability . Required for the localization of the long non-coding Xist RNA on the inactive chromosome X (Xi) and the subsequent initiation and maintenance of X-linked transcriptional gene silencing during X-inactivation . Plays a role as a RNA polymerase II (Pol II) holoenzyme transcription regulator .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Mouse monoclonal antibody to hnRNP U
Specificity	Recognizes endogenous levels of hnRNP U protein.
Antibody Type	Primary antibody
Immunogen	Purified recombinant human hnRNP U protein fragments expressed in E.coli.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 91 kD; Observed: 110 kD
Form/Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Alternative Names	HNRPU; SAFA; U21.1; Heterogeneous nuclear ribonucleoprotein U; hnRNP U; Scaffold attachment factor A; SAF-A; p120; pp120
Gene Symbol	HNRNPU
Entrez Gene	3192(Human)
SwissProt	Q00839(Human)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**hnRNP U Mouse Monoclonal Antibody(C3326)**

CAT. NO. AMA02938

DATA

Western blot analysis of hnRNP U expression in K562 (A), Jurkat (B), HeLa (C) whole cell lysates. (Predicted band size: 91 kD; Observed band size: 110 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.