

DATASHEET

Beclin-1 Mouse Monoclonal Antibody(C3338)

CAT. NO. AMA02950

KEY FEATURES

Target	Beclin-1	Source / Host	Mouse
Reactivity	Human, Mouse, Rat, Monkey, Hamster	Clonality	Monoclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.	Storage	at-20°C

BACKGROUND

Plays a central role in autophagy . Acts as a core subunit of the PI3K complex that mediates formation of phosphatidylinositol 3-phosphate; different complex forms are believed to play a role in multiple membrane trafficking pathways: PI3KC3-C1 is involved in initiation of autophagosomes and PI3KC3-C2 in maturation of autophagosomes and endocytosis. Involved in regulation of degradative endocytic trafficking and required for the abscission step in cytokinesis, probably in the context of PI3KC3-C2 . Essential for the formation of PI3KC3-C2 but not PI3KC3-C1 PI3K complex forms. Involved in endocytosis . May play a role in antiviral host defense.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Mouse monoclonal antibody to Beclin-1
Specificity	Recognizes endogenous levels of Beclin-1 protein.
Antibody Type	Primary antibody
Immunogen	Purified recombinant human Beclin 1 protein fragments expressed in E.coli.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 52 kD; Observed: 55 kD
Form/Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Alternative Names	GT197; Beclin-1; Coiled-coil myosin-like BCL2-interacting protein; Protein GT197
Gene Symbol	BECN1
Entrez Gene	8678(Human); 56208(Mouse); 114558(Rat)
SwissProt	Q14457(Human); O88597(Mouse); Q91XJ1(Rat)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**Beclin-1 Mouse Monoclonal Antibody(C3338)**

CAT. NO. AMA02950

DATA

Western blot analysis of Beclin-1 expression in HeLa (A), NIH3T3 (B), C6 (C), COS7 (D), CHOK1 (E), rat brain (F), Jurkat (G) whole cell lysates. (Predicted band size: 52 kD; Observed band size: 55 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.