

## DATASHEET

# RCC1 Mouse Monoclonal Antibody(C3714)

CAT. NO. AMA03326

### KEY FEATURES

Target	RCC1	Source / Host	Mouse
Reactivity	Human, Mouse	Clonality	Monoclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.	Storage	at-20°C

### BACKGROUND

Guanine-nucleotide releasing factor that promotes the exchange of Ran-bound GDP by GTP, and thereby plays an important role in RAN-mediated functions in nuclear import and mitosis . Contributes to the generation of high levels of chromosome-associated, GTP-bound RAN, which is important for mitotic spindle assembly and normal progress through mitosis . Via its role in maintaining high levels of GTP-bound RAN in the nucleus, contributes to the release of cargo proteins from importins after nuclear import . Involved in the regulation of onset of chromosome condensation in the S phase . Binds both to the nucleosomes and double-stranded DNA .

### APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:1000 - 1:2000
----	-----------------

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

### PRODUCT OVERVIEW


Description	Mouse monoclonal antibody to RCC1
Specificity	Recognizes endogenous levels of RCC1 protein.
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human RCC1. The exact sequence is proprietary.
Purification	This antibody is purified through a protein G column.
Molecular Weight	Predicted: 44 kD; Observed: 44 kD
Form/Buffer	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	CHC1; Regulator of chromosome condensation; Cell cycle regulatory protein; Chromosome condensation protein 1
Gene Symbol	RCC1
Entrez Gene	1104(Human); 100088(Mouse)
SwissProt	P18754(Human); Q8VE37(Mouse)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****RCC1 Mouse Monoclonal Antibody(C3714)**

CAT. NO. AMA03326

**DATA**

Western blot analysis of RCC1 expression in HeLa (A), HL60 (B), COS7 (C), C2C12 (D) whole cell lysates. (Predicted band size: 44 kD; Observed band size: 44 kD)

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.