

DATASHEET

CD73 Mouse Monoclonal Antibody(C3727)

CAT. NO. AMA03339

KEY FEATURES

Target	CD73	Source / Host	Mouse
Reactivity	Human	Clonality	Monoclonal
Applications	WB, IHC	Conjugation	Unconjugated
Form / Buffer	Mouse IgG1. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.	Storage	at-20°C

BACKGROUND

Catalyzes the hydrolysis of nucleotide monophosphates, releasing inorganic phosphate and the corresponding nucleoside, with AMP being the preferred substrate. Shows a preference for ribonucleotide monophosphates over their equivalent deoxyribose forms. Other substrates include IMP, UMP, GMP, CMP, dAMP, dCMP, dTMP, NAD and NMN.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:1000 - 1:2000
IHC	1:50 - 1:200

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

Description	Mouse monoclonal antibody to CD73
Specificity	Recognizes endogenous levels of CD73 protein.
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human CD73. The exact sequence is proprietary.
Purification	This antibody is purified through a protein G column.
Molecular Weight	Predicted: 63 kD; Observed: 70 kD
Form/Buffer	Mouse IgG1. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	NT5; NTE; 5'-nucleotidase; 5'-NT; Ecto-5'-nucleotidase; CD73
Gene Symbol	NT5E
Entrez Gene	4907(Human)
SwissProt	P21589(Human)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.


*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET


CD73 Mouse Monoclonal Antibody(C3727)

CAT. NO. AMA03339

DATA



Western blot analysis of CD73 expression in MDAMB231 (A) whole cell lysates. (Predicted band size: 63 kD; Observed band size: 70 kD)



Immunohistochemical analysis of CD73 staining in human kidney formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.