

DATASHEET

PAX6 Mouse Monoclonal Antibody(C3746)

CAT. NO. AMA03358

KEY FEATURES

Target	PAX6	Source / Host	Mouse
Reactivity	Human	Clonality	Monoclonal
Applications	WB, IF/ICC	Conjugation	Unconjugated
Form / Buffer	Mouse IgM kappa. Supplied in crude ascites with 0.01% sodium azide.	Storage	at-20°C

BACKGROUND

Transcription factor with important functions in the development of the eye, nose, central nervous system and pancreas. Required for the differentiation of pancreatic islet alpha cells . Competes with PAX4 in binding to a common element in the glucagon, insulin and somatostatin promoters. Regulates specification of the ventral neuron subtypes by establishing the correct progenitor domains . Acts as a transcriptional repressor of NFATC1-mediated gene expression .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:4000
IF/ICC	1:10 - 1:50

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Mouse monoclonal antibody to PAX6
Specificity	Recognizes endogenous levels of PAX6 protein.
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human PAX6. The exact sequence is proprietary.
Molecular Weight	Predicted: 46 kD; Observed: 50 kD
Form/Buffer	Mouse IgM kappa. Supplied in crude ascites with 0.01% sodium azide.
Alternative Names	AN2; Paired box protein Pax-6; Aniridia type II protein; Oculorhombin
Gene Symbol	PAX6
Entrez Gene	5080(Human)
SwissProt	P26367(Human)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.


*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**PAX6 Mouse Monoclonal Antibody(C3746)**

CAT. NO. AMA03358

DATA

Western blot analysis of PAX6 expression in 293 (A) whole cell lysates. (Predicted band size: 46 kD; Observed band size: 50 kD)



Immunofluorescent analysis of PAX6 staining in HeLa cells. Formalin-fixed cells were permeabilized with 0.1% Triton X-100 in TBS for 5-10 minutes and blocked with 3% BSA-PBS for 30 minutes at room temperature. Cells were probed with the primary antibody in 3% BSA-PBS and incubated overnight at 4 °C in a humidified chamber. Cells were washed with PBST and incubated with a AREX® Fluor 488 -conjugated secondary antibody (green) in PBS at room temperature in the dark. Phalloidin - AREX® Fluor 555 was used to stain Actin filaments (red). DAPI was used to stain the cell nuclei (blue).

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.