

## DATASHEET

# Aconitase 2 Mouse Monoclonal Antibody(C3979)

CAT. NO. AMA03591

### KEY FEATURES

Target	Aconitase 2	Source / Host	Mouse
Reactivity	Human	Clonality	Monoclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Mouse IgG1. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.	Storage	at-20°C

### BACKGROUND

Catalyzes the isomerization of citrate to isocitrate via cis-aconitate.

### APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
----	----------------

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

### PRODUCT OVERVIEW


Description	Mouse monoclonal antibody to Aconitase 2
Specificity	Recognizes endogenous levels of Aconitase 2 protein.
Antibody Type	Primary antibody
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Aconitase 2. The exact sequence is proprietary.
Purification	This antibody is purified through a protein G column.
Molecular Weight	Predicted: 85 kD; Observed: 85 kD
Form/Buffer	Mouse IgG1. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	Aconitate hydratase, mitochondrial; Aconitase; Citrate hydro-lyase
Gene Symbol	ACO2
Entrez Gene	50(Human)
SwissProt	Q99798(Human)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****Aconitase 2 Mouse Monoclonal Antibody(C3979)**

CAT. NO. AMA03591

**DATA**

Western blot analysis of Aconitase 2 expression in human heart (A) whole cell lysates.  
(Predicted band size: 85 kD; Observed band size: 85 kD)

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.