

## DATASHEET

# CD96 Rat Monoclonal Antibody(3.3)

CAT. NO. AMA03705

### KEY FEATURES

Target	CD96	Source / Host	Rat
Reactivity	Mouse	Clonality	Monoclonal
Applications	IF/ICC, FC	Conjugation	Unconjugated
Form / Buffer	Rat IgG1 kappa. Liquid in PBS, pH 7.3, and 0.02% sodium azide.		Storage at-20°C

### BACKGROUND

The protein encoded by this gene belongs to the immunoglobulin superfamily. It is a type I membrane protein. The protein may play a role in the adhesive interactions of activated T and NK cells during the late phase of the immune response. It may also function in antigen presentation. Alternative splicing generates multiple transcript variants encoding distinct isoforms.

### APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

IF/ICC	1:50 - 1:200
FC	1:500 - 1:2000

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

### PRODUCT OVERVIEW

Description	Rat monoclonal antibody to CD96
Specificity	Recognizes mouse CD96
Antibody Type	Primary antibody
Immunogen	Proprietary immunogen
Purification	The antibody was purified by affinity chromatography.
Form/Buffer	Rat IgG1 kappa. Liquid in PBS, pH 7.3, and 0.02% sodium azide.
Alternative Names	T-cell surface protein tactile; Cell surface antigen CD96; T cell-activated increased late expression protein; CD96
Gene Symbol	CD96
Entrez Gene	84544(Mouse)
SwissProt	Q3U0X8(Mouse)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****CD96 Rat Monoclonal Antibody(3.3)**

CAT. NO. AMA03705

**| DATA****| STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**| NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.