

DATASHEET

CD155 Rat Monoclonal Antibody(4.24.1)

CAT. NO. AMA03716

KEY FEATURES

Target	CD155	Source / Host	Rat
Reactivity	Mouse	Clonality	Monoclonal
Applications	IF/ICC, FC	Conjugation	Unconjugated
Form / Buffer	Rat IgG2a kappa. Liquid in PBS, pH 7.3, and 0.02% sodium azide.	Storage	at-20°C

BACKGROUND

The protein encoded by this gene is a transmembrane glycoprotein belonging to the immunoglobulin superfamily. The external domain mediates cell attachment to the extracellular matrix molecule vitronectin, while its intracellular domain interacts with the dynein light chain Tctex-1/DYNLT1. The gene is specific to the primate lineage, and serves as a cellular receptor for poliovirus in the first step of poliovirus replication. Multiple transcript variants encoding different isoforms have been found for this gene.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

PRODUCT OVERVIEW

Description	Rat monoclonal antibody to CD155
Specificity	Recognizes mouse CD155
Antibody Type	Primary antibody
Immunogen	CD155 transfectants
Purification	The antibody was purified by affinity chromatography.
Form/Buffer	Rat IgG2a kappa. Liquid in PBS, pH 7.3, and 0.02% sodium azide.
Alternative Names	Cd155; Poliovirus receptor
Gene Symbol	PVR
Entrez Gene	52118(Mouse)
SwissProt	Q91WP1(Mouse)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET

CD155 Rat Monoclonal Antibody(4.24.1)

CAT. NO. AMA03716

| DATA

| STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

| NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.