

DATASHEET

TCR V alpha 24-J alpha 18 Mouse Monoclonal Antibody(6B11)

CAT. NO. AMA03730

KEY FEATURES

Target	TCR V alpha 24-J alpha 18	Source / Host	Mouse
Reactivity	Human	Clonality	Monoclonal
Applications	IF/ICC, FC	Conjugation	Unconjugated
Form / Buffer	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, and 0.02% sodium azide.	Storage	at-20°C

BACKGROUND

This gene encodes a tumor suppressor protein containing transcriptional activation, DNA binding, and oligomerization domains. The encoded protein responds to diverse cellular stresses to regulate expression of target genes, thereby inducing cell cycle arrest, apoptosis, senescence, DNA repair, or changes in metabolism. Mutations in this gene are associated with a variety of human cancers, including hereditary cancers such as Li-Fraumeni syndrome. Alternative splicing of this gene and the use of alternate promoters result in multiple transcript variants and isoforms. Additional isoforms have also been shown to result from the use of alternate translation initiation codons from identical transcript variants (PMIDs: 12032546, 20937277).

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

IF/ICC	1:50 - 1:200
FC	1:500 - 1:2000

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

Description	Mouse monoclonal antibody to TCR V alpha 24-J alpha 18
Specificity	Recognizes human TCR V alpha 24-J alpha 18
Antibody Type	Primary antibody
Immunogen	Peptide epitope in the CDR3 region of V alpha 24
Purification	The antibody was purified by affinity chromatography.
Form/Buffer	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, and 0.02% sodium azide.
Alternative Names	iNKT cell; iNKT; invariant natural killer T cel; invariant NK T cell; invariant NKT cell
Gene Symbol	TRA
Entrez Gene	6955(Human)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**TCR V alpha 24-J alpha 18 Mouse Monoclonal Antibody(6B11)****CAT. NO. AMA03730****| DATA****| STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

| NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.