

DATASHEET

CD183 Armenian Hamster Monoclonal Antibody(CXCR3-173)

CAT. NO. AMA03769

KEY FEATURES

Target	CD183	Source / Host	Armenian Hamster
Reactivity	Mouse	Clonality	Monoclonal
Applications	IF/ICC, FC	Conjugation	Unconjugated
Form / Buffer	Armenian Hamster IgG. Liquid in PBS, Storage pH 7.3, and 0.02% sodium azide.		at-20°C

BACKGROUND

This gene encodes a G protein-coupled receptor with selectivity for three chemokines, termed CXCL9/Mig (monokine induced by interferon-g), CXCL10/IP10 (interferon-g-inducible 10 kDa protein) and CXCL11/I-TAC (interferon-inducible T cell a-chemoattractant). Binding of chemokines to this protein induces cellular responses that are involved in leukocyte traffic, most notably integrin activation, cytoskeletal changes and chemotactic migration. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. One of the isoforms (CXCR3-B) shows high affinity binding to chemokine, CXCL4/PF4 (PMID:12782716).

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

IF/ICC	1:50 - 1:200
FC	1:500 - 1:2000

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

Description	Armenian Hamster monoclonal antibody to CD183
Specificity	Recognizes mouse CD183
Antibody Type	Primary antibody
Immunogen	Mouse N-terminus of CXCR3
Purification	The antibody was purified by affinity chromatography.
Form/Buffer	Armenian Hamster IgG. Liquid in PBS, pH 7.3, and 0.02% sodium azide.
Alternative Names	GPR9; C-X-C chemokine receptor type 3; CXC-R3; CXCR-3; CKR-L2; G protein-coupled receptor 9; Interferon-inducible protein 10 receptor; IP-10 receptor; CD183
Gene Symbol	CXCR3
Entrez Gene	12766(Mouse)
SwissProt	O88410(Mouse)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**CD183 Armenian Hamster Monoclonal Antibody(CXCR3-173)**

CAT. NO. AMA03769

| DATA**| STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

| NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.