

DATASHEET

CD11a Rat Monoclonal Antibody(I21/7)

CAT. NO. AMA03807

KEY FEATURES

Target	CD11a	Source / Host	Rat
Reactivity	Mouse	Clonality	Monoclonal
Applications	IF/ICC, FC	Conjugation	Unconjugated
Form / Buffer	Rat IgG2a kappa. Liquid in PBS, pH 7.3, and 0.02% sodium azide.	Storage	at-20°C

BACKGROUND

ITGAL encodes the integrin alpha L chain. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. This I-domain containing alpha integrin combines with the beta 2 chain (ITGB2) to form the integrin lymphocyte function-associated antigen-1 (LFA-1), which is expressed on all leukocytes. LFA-1 plays a central role in leukocyte intercellular adhesion through interactions with its ligands, ICAMs 1-3 (intercellular adhesion molecules 1 through 3), and also functions in lymphocyte costimulatory signaling. Two transcript variants encoding different isoforms have been found for this gene.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

PRODUCT OVERVIEW

Description	Rat monoclonal antibody to CD11a
Specificity	Recognizes mouse CD11a
Antibody Type	Primary antibody
Immunogen	C57BL/6 mouse splenic secondary cytotoxic T cells
Purification	The antibody was purified by affinity chromatography.
Form/Buffer	Rat IgG2a kappa. Liquid in PBS, pH 7.3, and 0.02% sodium azide.
Alternative Names	CD11A; Integrin alpha-L; CD11 antigen-like family member A; Leukocyte adhesion glycoprotein LFA-1 alpha chain; LFA-1A; Leukocyte function-associated molecule 1 alpha chain; CD11a
Gene Symbol	ITGAL
Entrez Gene	16408(Mouse)
SwissProt	P24063(Mouse)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**CD11a Rat Monoclonal Antibody(I21/7)**

CAT. NO. AMA03807

| DATA**| STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

| NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.