

DATASHEET

IL-4 Rat Monoclonal Antibody(MP4-25D2)

CAT. NO. AMA03832

KEY FEATURES

Target	IL-4	Source / Host	Rat
Reactivity	Human	Clonality	Monoclonal
Applications	IF/ICC, FC	Conjugation	Unconjugated
Form / Buffer	Rat IgG1 kappa. Liquid in PBS, pH 7.3, and 0.02% sodium azide.		Storage at-20°C

BACKGROUND

Cytokine secreted primarily by mast cells, T-cells, eosinophils, and basophils that plays a role in regulating antibody production, hematopoiesis and inflammation, and the development of effector T-cell responses . Induces the expression of class II MHC molecules on resting B-cells. Enhances both secretion and cell surface expression of IgE and IgG1 . Also regulates the expression of the low affinity Fc receptor for IgE (CD23) on both lymphocytes and monocytes . Positively regulates IL31RA expression in macrophages. Stimulates autophagy in dendritic cells by interfering with mTORC1 signaling and through the induction of RUFY4. In addition, plays a critical role in higher functions of the normal brain, such as memory and learning .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

IF/ICC	1:50 - 1:200
FC	1:500 - 1:2000

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

Description	Rat monoclonal antibody to IL-4
Specificity	Recognizes human IL-4
Antibody Type	Primary antibody
Immunogen	CHO-expressed recombinant human IL-4
Purification	The antibody was purified by affinity chromatography.
Form/Buffer	Rat IgG1 kappa. Liquid in PBS, pH 7.3, and 0.02% sodium azide.
Alternative Names	Interleukin-4; IL-4; B-cell stimulatory factor 1; BSF-1; Binetrakin; Lymphocyte stimulatory factor 1; Pitrakinra
Gene Symbol	IL4
Entrez Gene	3565(Human)
SwissProt	P05112(Human)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET

IL-4 Rat Monoclonal Antibody(MP4-25D2)

CAT. NO. AMA03832

| DATA

| STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

| NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.