

## DATASHEET

# TNF alpha Armenian Hamster Monoclonal Antibody(TN3-19.12)

CAT. NO. AMA03861

### KEY FEATURES

Target	TNF alpha	Source / Host	Armenian Hamster
Reactivity	Mouse	Clonality	Monoclonal
Applications	IF/ICC, FC	Conjugation	Unconjugated
Form / Buffer	Armenian Hamster IgG. Liquid in PBS, Storage pH 7.3, and 0.02% sodium azide.		at-20°C

### BACKGROUND

This gene encodes a multifunctional proinflammatory cytokine that belongs to the tumor necrosis factor (TNF) superfamily. This cytokine is mainly secreted by macrophages. It can bind to, and thus functions through its receptors TNFRSF1A/TNFR1 and TNFRSF1B/TNFR2. This cytokine is involved in the regulation of a wide spectrum of biological processes including cell proliferation, differentiation, apoptosis, lipid metabolism, and coagulation. This cytokine has been implicated in a variety of diseases, including autoimmune diseases, insulin resistance, psoriasis, rheumatoid arthritis ankylosing spondylitis, tuberculosis, autosomal dominant polycystic kidney disease, and cancer. Mutations in this gene affect susceptibility to cerebral malaria, septic shock, and Alzheimer disease.

### APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

IF/ICC	1:50 - 1:200
FC	1:500 - 1:2000

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

### PRODUCT OVERVIEW

Description	Armenian Hamster monoclonal antibody to TNF alpha
Specificity	Recognizes mouse TNF alpha
Antibody Type	Primary antibody
Immunogen	E. coli-expressed, recombinant mouse TNF-α
Purification	The antibody was purified by affinity chromatography.
Form/Buffer	Armenian Hamster IgG. Liquid in PBS, pH 7.3, and 0.02% sodium azide.
Alternative Names	TNFA; TNFSF2; Tumor necrosis factor; Cachectin; TNF-alpha; Tumor necrosis factor ligand superfamily member 2; TNF-a
Gene Symbol	TNF
Entrez Gene	21926(Mouse)
SwissProt	P06804(Mouse)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****TNF alpha Armenian Hamster Monoclonal Antibody(TN3-19.12)**

CAT. NO. AMA03861

**| DATA****| STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**| NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.