

## DATASHEET

# V gamma 2 TCR Armenian Hamster Monoclonal Antibody(UC3-10A6)

CAT. NO. AMA03869

### KEY FEATURES

Target	V gamma 2 TCR	Source / Host	Armenian Hamster
Reactivity	Mouse	Clonality	Monoclonal
Applications	IF/ICC, FC	Conjugation	Unconjugated
Form / Buffer	Armenian Hamster IgG. Liquid in PBS, Storage pH 7.3, and 0.02% sodium azide.		at-20°C

### BACKGROUND

This gene encodes a tumor suppressor protein containing transcriptional activation, DNA binding, and oligomerization domains. The encoded protein responds to diverse cellular stresses to regulate expression of target genes, thereby inducing cell cycle arrest, apoptosis, senescence, DNA repair, or changes in metabolism. Mutations in this gene are associated with a variety of human cancers, including hereditary cancers such as Li-Fraumeni syndrome. Alternative splicing of this gene and the use of alternate promoters result in multiple transcript variants and isoforms. Additional isoforms have also been shown to result from the use of alternate translation initiation codons from identical transcript variants (PMIDs: 12032546, 20937277).

### APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

IF/ICC	1:50 - 1:200
FC	1:500 - 1:2000

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

### PRODUCT OVERVIEW

Description	Armenian Hamster monoclonal antibody to V gamma 2 TCR
Specificity	Recognizes mouse V gamma 2 TCR
Antibody Type	Primary antibody
Immunogen	G8 mouse T cells
Purification	The antibody was purified by affinity chromatography.
Form/Buffer	Armenian Hamster IgG. Liquid in PBS, pH 7.3, and 0.02% sodium azide.
Alternative Names	V gamma 2 TCR; Vy2 TCR; Trgv2
Gene Symbol	Tcrg
Entrez Gene	110067(Mouse)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****V gamma 2 TCR Armenian Hamster Monoclonal Antibody(UC3-10A6)****CAT. NO. AMA03869****| DATA****| STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**| NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.