

DATASHEET

Cleaved Caspase 2 Rabbit Polyclonal Antibody

CAT. NO. APA06729

KEY FEATURES

Target	Cleaved Caspase 2	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

Is a regulator of the cascade of caspases responsible for apoptosis execution . Might function by either activating some proteins required for cell death or inactivating proteins necessary for cell survival . Associates with PIDD1 and CRADD to form the PIDDosome, a complex that activates CASP2 and triggers apoptosis in response to genotoxic stress .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to Cleaved Caspase 2
Specificity	Recognizes endogenous levels of Cleaved Caspase 2 protein.
Antibody Type	Primary antibody
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Cleaved Caspase 2. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 50 kD; Observed: 34 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	ICH1; NEDD2; Caspase-2; CASP-2; Neural precursor cell expressed developmentally down-regulated protein 2; NEDD-2; Protease ICH-1
Gene Symbol	CASP2
Entrez Gene	835(Human); 12366(Mouse); 64314(Rat)
SwissProt	P42575(Human); P29594(Mouse); P55215(Rat)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**Cleaved Caspase 2 Rabbit Polyclonal Antibody**

CAT. NO. APA06729

DATA

Western blot analysis of Cleaved Caspase 2 expression in HepG2 (A), mouse liver (B), mouse brain (C), rat liver (D), rat brain (E), rat heart (F) whole cell lysates. (Predicted band size: 50 kD; Observed band size: 34 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.