

DATASHEET

GRLF1 (Phospho-Y1087) Rabbit Polyclonal Antibody

CAT. NO. APA07051

KEY FEATURES

Target	GRLF1 (Phospho-Y1087)	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat, Dog, Pig	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

Rho GTPase-activating protein (GAP) . Binds several acidic phospholipids which inhibits the Rho GAP activity to promote the Rac GAP activity . This binding is inhibited by phosphorylation by PRKCA . Involved in cell differentiation as well as cell adhesion and migration, plays an important role in retinal tissue morphogenesis, neural tube fusion, midline fusion of the cerebral hemispheres and mammary gland branching morphogenesis . Transduces signals from p21-ras to the nucleus, acting via the ras GTPase-activating protein (GAP) . Transduces SRC-dependent signals from cell-surface adhesion molecules, such as laminin, to promote neurite outgrowth. Regulates axon outgrowth, guidance and fasciculation .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to GRLF1 (Phospho-Y1087)
Specificity	Recognizes endogenous levels of GRLF1 protein only when phosphorylated at Y1087.
Antibody Type	Primary antibody
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y1087 of human GRLF1 protein. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 170 kD; Observed: 190 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	GRF1; GRLF1; KIAA1722; Rho GTPase-activating protein 35; Glucocorticoid receptor DNA-binding factor 1; Glucocorticoid receptor repression factor 1; GRF-1; Rho GAP p190A; p190-A
Gene Symbol	ARHGAP35
Entrez Gene	2909(Human); 232906(Mouse)
SwissProt	Q9NRY4(Human); Q91YM2(Mouse); P81128(Rat)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**GRLF1 (Phospho-Y1087) Rabbit Polyclonal Antibody**

CAT. NO. APA07051

DATA

Western blot analysis of GRLF1 (Phospho-Y1087) expression in mouse testis (A), MCF7 (B), HUT78 (C), K562 (D) whole cell lysates. (Predicted band size: 170 kD; Observed band size: 190 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.