

DATASHEET

HEXIM2 Rabbit Polyclonal Antibody

CAT. NO. APA08275

KEY FEATURES

| | | | |
|---------------|---|---------------|---------------------|
| Target | HEXIM2 | Source / Host | Rabbit |
| Reactivity | Human, Rat, Bovine | Clonality | Polyclonal |
| Applications | WB | Conjugation | Unconjugated |
| Form / Buffer | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. | | Storage at -20°C |

BACKGROUND

Transcriptional regulator which functions as a general RNA polymerase II transcription inhibitor . Core component of the 7SK RNP complex: in cooperation with 7SK snRNA sequesters P-TEFb in a large inactive 7SK snRNP complex preventing RNA polymerase II phosphorylation and subsequent transcriptional elongation .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

| | |
|----|----------------|
| WB | 1:500 - 1:1000 |
|----|----------------|

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


| | |
|-------------------|---|
| Description | Rabbit polyclonal antibody to HEXIM2 |
| Specificity | Recognizes endogenous levels of HEXIM2 protein. |
| Antibody Type | Primary antibody |
| Immunogen | KLH-conjugated synthetic peptide encompassing a sequence within the center region of human HEXIM2. The exact sequence is proprietary. |
| Purification | The antibody was purified by immunogen affinity chromatography. |
| Molecular Weight | Predicted: 32 kD; Observed: 42 kD |
| Form/Buffer | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide. |
| Alternative Names | Protein HEXIM2; Hexamethylene bis-acetamide-inducible protein 2 |
| Gene Symbol | HEXIM2 |
| Entrez Gene | 124790(Human) |
| SwissProt | Q96MH2(Human) |

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**HEXIM2 Rabbit Polyclonal Antibody**

CAT. NO. APA08275

DATA

Western blot analysis of HEXIM2 expression in HEK293T (A), rat kidney (B) whole cell lysates. (Predicted band size: 32 kD; Observed band size: 42 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.