

**DATASHEET**

**ABCF3 Rabbit Polyclonal Antibody**

CAT. NO. APA08477

**KEY FEATURES**

Target	ABCF3	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat, Bovine, Pig	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

**BACKGROUND**

Displays an antiviral effect against flaviviruses such as west Nile virus (WNV) in the presence of OAS1B.

**APPLICATION**

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
----	----------------

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

**PRODUCT OVERVIEW**


Description	Rabbit polyclonal antibody to ABCF3
Specificity	Recognizes endogenous levels of ABCF3 protein.
Antibody Type	Primary antibody
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human ABCF3. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 78; Observed: 80 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	ATP-binding cassette sub-family F member 3
Gene Symbol	ABCF3
Entrez Gene	55324(Human); 27406(Mouse); 287982(Rat)
SwissProt	Q9NUQ8(Human); Q8K268(Mouse); Q66H39(Rat)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****ABCF3 Rabbit Polyclonal Antibody**

CAT. NO. APA08477

**DATA**

Western blot analysis of ABCF3 expression in A549 (A), H9C2 (B) whole cell lysates. (Predicted band size: 78; 79 kD; Observed band size: 80 kD)

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.