

## DATASHEET

# ABHD7 Rabbit Polyclonal Antibody

CAT. NO. APA09671

### KEY FEATURES

Target	ABHD7	Source / Host	Rabbit
Reactivity	Human	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

### BACKGROUND

Epoxide hydrolase 4 (EPHX4), also known as Abhydrolase domain-containing protein 7, Epoxide hydrolase-related protein. Belongs to the AB hydrolase superfamily. Epoxide hydrolase family Subcellular location: Membrane

### APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
----	----------------

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

### PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to ABHD7
Specificity	Recognizes endogenous levels of ABHD7 protein.
Antibody Type	Primary antibody
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human ABHD7. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 42 kD; Observed: 42 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	ABHD7; EPHXRP; Epoxide hydrolase 4; Abhydrolase domain-containing protein 7; Epoxide hydrolase-related protein
Gene Symbol	EPHX4
Entrez Gene	253152(Human)
SwissProt	Q8IUS5(Human)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****ABHD7 Rabbit Polyclonal Antibody**

CAT. NO. APA09671

**DATA**

Western blot analysis of ABHD7 expression in HEK293T (A), HepG2 (B) whole cell lysates.  
(Predicted band size: 42 kD; Observed band size: 42 kD)

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.