

DATASHEET

cTnI (Phospho-S22/S23) Rabbit Polyclonal Antibody

CAT. NO. APA09793

KEY FEATURES

Target	cTnI (Phospho-S22/S23)	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat, Pig, Rabbit	Clonality	Polyclonal
Applications	WB, IHC	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.	Storage	at-20°C

BACKGROUND

Troponin I is the inhibitory subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
IHC	1:100 - 1:200

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

Description	Rabbit polyclonal antibody to cTnI (Phospho-S22/S23)
Specificity	Recognizes endogenous levels of cTnI protein only when phosphorylated at S22/S23.
Antibody Type	Primary antibody
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S22/S23 of human cTnI protein. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 24 kD; Observed: 28 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	TNNI3; TNNC1; Troponin I, cardiac muscle; Cardiac troponin I
Gene Symbol	TNNI3
Entrez Gene	7137(Human); 21954(Mouse); 29248(Rat)
SwissProt	P19429(Human); P48787(Mouse); P23693(Rat)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.


*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET


cTnI (Phospho-S22/S23) Rabbit Polyclonal Antibody

CAT. NO. APA09793

DATA



Western blot analysis of cTnI (Phospho-S22/S23) expression in mouse heart (A), mouse kidney (B), rat heart (C), rat kidney (D) whole cell lysates. (Predicted band size: 24 kD; Observed band size: 28 kD)



Immunohistochemical analysis of cTnI (Phospho-S22/S23) staining in human heart formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.