

DATASHEET

FAM84B Rabbit Polyclonal Antibody

CAT. NO. APA10190

KEY FEATURES

Target	FAM84B	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat, Bovine, Monkey	Clonality	Polyclonal
Applications	WB, IHC	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

Protein LRATD2 (LRATD2), also known as Breast cancer membrane protein 101, LRAT domain-containing 2, Protein FAM84B. Belongs to the LRATD family Tissue specificity: Expressed in esophageal squamous cell carcinomas

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
IHC	1:100 - 1:200

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

Description	Rabbit polyclonal antibody to FAM84B
Specificity	Recognizes endogenous levels of FAM84B protein.
Antibody Type	Primary antibody
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human FAM84B. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 34 kD; Observed: 36 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	BCMP101; NSE2; Protein FAM84B; Breast cancer membrane protein 101; Protein NSE2
Gene Symbol	FAM84B
Entrez Gene	157638(Human)
SwissProt	Q96KN1(Human)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.


*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET


FAM84B Rabbit Polyclonal Antibody

CAT. NO. APA10190

DATA



Western blot analysis of FAM84B expression in HUVEC (A) whole cell lysates. (Predicted band size: 34 kD; Observed band size: 36 kD)



Immunohistochemical analysis of FAM84B staining in human breast cancer formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.