

DATASHEET

Tachykinin Receptor 3 Rabbit Polyclonal Antibody

CAT. NO. APA12312

KEY FEATURES

Target	Tachykinin Receptor 3	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

Receptor for the tachykinin neuromedin-K (neurokinin B), also able to bind and respond to tachynins substance K/neurokinin A and substance P , also able to bind and respond to tachynins substance K/neurokinin A and substance P . The rank order of affinity of this receptor to tachykinins is: neuromedin-K > substance K and substance P . Neuromedin-K binding to its receptor triggers G protein-coupled receptor signaling via activation of G(q) and phosphatidylinositol hydrolysis by phospholipase C . Neuromedin-K binding also triggers signaling via activation of adenylate cyclase activity which results in increased intracellular levels of cyclic AMP (cAMP) .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to Tachykinin Receptor 3
Specificity	Recognizes endogenous levels of Tachykinin Receptor 3 protein.
Antibody Type	Primary antibody
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Tachykinin Receptor 3. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 52 kD; Observed: 52 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	NK3R; TAC3R; Neuromedin-K receptor; NKR; NK-3 receptor; NK-3R; Neurokinin B receptor; Tachykinin receptor 3
Gene Symbol	TACR3
Entrez Gene	6870(Human); 21338(Mouse); 24808(Rat)
SwissProt	P29371(Human); P47937(Mouse); P16177(Rat)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**Tachykinin Receptor 3 Rabbit Polyclonal Antibody**

CAT. NO. APA12312

DATA

Western blot analysis of Tachykinin Receptor 3 expression in NIH3T3 (A) whole cell lysates. (Predicted band size: 52 kD; Observed band size: 52 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.