

**DATASHEET**

**FITC-conjugated Ovalbumin Rabbit Polyclonal Antibody**

CAT. NO. APA12509

**KEY FEATURES**

Target	Ovalbumin	Source / Host	Rabbit
Reactivity	Chicken	Clonality	Polyclonal
Applications	ELISA	Conjugation	FITC
Form / Buffer	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.	Storage	at-20°C

**BACKGROUND**

Ovalbumin (OVA) is the major protein of egg white (~45 kDa) and a model antigen extensively used in immunology research, particularly in murine models of allergy, asthma, and immune tolerance. Anti-ovalbumin antibodies are used in ELISA, WB, and immunoassay for tracking OVA in experimental systems and detecting egg protein contamination in food.

**APPLICATION**

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

ELISA	Use at an assay dependent dilution.
-------	-------------------------------------

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

**PRODUCT OVERVIEW**

Description	Rabbit polyclonal antibody to Ovalbumin-FITC labeled
Specificity	Reacts only with Ovalbumin.
Antibody Type	Primary antibody
Immunogen	Full length Ovalbumin
Purification	Protein G purified
Form/Buffer	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET**

**FITC-conjugated Ovalbumin Rabbit Polyclonal Antibody**

**CAT. NO. APA12509**

**| DATA**

**| STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**| NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.