

DATASHEET

Biotin-conjugated Human Serum Albumin Rabbit Polyclonal Antibody

CAT. NO. APA12516

KEY FEATURES

Target	Human Serum Albumin	Source / Host	Rabbit
Reactivity	Human	Clonality	Polyclonal
Applications	ELISA	Conjugation	Biotin
Form / Buffer	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.	Storage	at-20°C

BACKGROUND

Binds water, Ca(2+), Na(+), K(+), fatty acids, hormones, bilirubin and drugs (Probable). Its main function is the regulation of the colloidal osmotic pressure of blood (Probable). Major zinc transporter in plasma, typically binds about 80% of all plasma zinc , Na(+), K(+), fatty acids, hormones, bilirubin and drugs (Probable). Its main function is the regulation of the colloidal osmotic pressure of blood (Probable). Major zinc transporter in plasma, typically binds about 80% of all plasma zinc . Major calcium and magnesium transporter in plasma, binds approximately 45% of circulating calcium and magnesium in plasma . Potentially has more than two calcium-binding sites and might additionally bind calcium in a non-specific manner .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

ELISA	Use at an assay dependent dilution.
-------	-------------------------------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

Description	Rabbit polyclonal antibody to Human Serum Albumin-Biotin labeled
Specificity	Reacts only with Human Serum Albumin.
Antibody Type	Primary antibody
Immunogen	Human Serum Albumin
Purification	Immunogen affinity purified
Form/Buffer	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.
Alternative Names	Serum albumin
Gene Symbol	ALB
Entrez Gene	213(Human)
SwissProt	P02768(Human)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**Biotin-conjugated Human Serum Albumin Rabbit Polyclonal Antibody****CAT. NO. APA12516****| DATA****| STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

| NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.