

DATASHEET

Fis1 Rabbit Polyclonal Antibody

CAT. NO. APA12806

KEY FEATURES

Target	Fis1	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB, IF/ICC, IP	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

Involved in the fragmentation of the mitochondrial network and its perinuclear clustering . Plays a minor role in the recruitment and association of the fission mediator dynamin-related protein 1 (DNM1L) to the mitochondrial surface and mitochondrial fission . May not be essential for the assembly of functional fission complexes and the subsequent membrane scission event . Also mediates peroxisomal fission . May act when the products of fission are directed toward mitochondrial homeostasis, mitophagy, or apoptosis . Can induce cytochrome c release from the mitochondrion to the cytosol, ultimately leading to apoptosis .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
IF/ICC	1:50 - 1:200
IP	1:10 - 1:50

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

Description	Rabbit polyclonal antibody to Fis1
Specificity	Recognizes endogenous levels of Fis1 protein.
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human Fis1
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 16 kD; Observed: 16 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	TTC11; Mitochondrial fission 1 protein; FIS1 homolog; hFis1; Tetratricopeptide repeat protein 11; TPR repeat protein 11
Gene Symbol	FIS1
Entrez Gene	51024(Human); 66437(Mouse); 288584(Rat)
SwissProt	Q9Y3D6(Human); Q9CQ92(Mouse); P84817(Rat)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.


*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET


Fis1 Rabbit Polyclonal Antibody

CAT. NO. APA12806

DATA



Western blot analysis of Fis1 expression in HepG2 (A), A549 (B), mouse spleen (C), mouse kidney (D) whole cell lysates. (Predicted band size: 16 kD; Observed band size: 16 kD)



Immunofluorescent analysis of Fis1 staining in MCF7 cells. Formalin-fixed cells were permeabilized with 0.1% Triton X-100 in TBS for 5-10 minutes and blocked with 3% BSA-PBS for 30 minutes at room temperature. Cells were probed with the primary antibody in 3% BSA-PBS and incubated overnight at 4 °C in a humidified chamber. Cells were washed with PBST and incubated with a DyLight 594-conjugated secondary antibody (red) in PBS at room temperature in the dark.

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.