

DATASHEET

COG2 Rabbit Polyclonal Antibody

CAT. NO. APA12838

KEY FEATURES

Target	COG2	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB, IHC	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

Required for normal Golgi morphology and function.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
IHC	1:50 - 1:200

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

Description	Rabbit polyclonal antibody to COG2
Specificity	Recognizes endogenous levels of COG2 protein.
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human COG2
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 83 kD; Observed: 88 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	LDLC; Conserved oligomeric Golgi complex subunit 2; COG complex subunit 2; Component of oligomeric Golgi complex 2; Low density lipoprotein receptor defect C-complementing protein
Gene Symbol	COG2
Entrez Gene	22796(Human); 76332(Mouse)
SwissProt	Q14746(Human); Q921L5(Mouse)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.


*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET


COG2 Rabbit Polyclonal Antibody

CAT. NO. APA12838

DATA



Western blot analysis of COG2 expression in HepG2 (A), MCF7 (B), mouse kidney (C), rat testis (D) whole cell lysates. (Predicted band size: 83 kD; Observed band size: 88 kD)



Immunohistochemical analysis of COG2 staining in rat brain formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.