

DATASHEET

SPR Rabbit Polyclonal Antibody

CAT. NO. APA13185

KEY FEATURES

Target	SPR	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

Catalyzes the final one or two reductions in tetra-hydrobiopterin biosynthesis to form 5,6,7,8-tetrahydrobiopterin.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to SPR
Specificity	Recognizes endogenous levels of SPR protein.
Antibody Type	Primary antibody
Immunogen	KLH-conjugated synthetic peptide of human SPR
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 28 kD; Observed: 31 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	Sepiapterin reductase; SPR
Gene Symbol	SPR
Entrez Gene	6697(Human); 29270(Rat)
SwissProt	P35270(Human); Q64105(Mouse); P18297(Rat)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**SPR Rabbit Polyclonal Antibody**

CAT. NO. APA13185

DATA

Western blot analysis of SPR expression in Hela (A), COS7 (B), MCF7 (C), mouse liver (D), rat kidney (E) whole cell lysates. (Predicted band size: 28 kD; Observed band size: 31 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.