

## DATASHEET

# STAR Rabbit Polyclonal Antibody

CAT. NO. APA13720

### KEY FEATURES

Target	STAR	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

### BACKGROUND

Plays a key role in steroid hormone synthesis by enhancing the metabolism of cholesterol into pregnenolone. Mediates the transfer of cholesterol from the outer mitochondrial membrane to the inner mitochondrial membrane where it is cleaved to pregnenolone.

### APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
----	----------------

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

### PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to STAR
Specificity	Recognizes endogenous levels of STAR protein.
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human STAR
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 31 kD; Observed: 31 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	STARD1; Steroidogenic acute regulatory protein mitochondrial; StAR; START domain-containing protein 1; StARD1
Gene Symbol	STAR
Entrez Gene	6770(Human); 20845(Mouse); 25557(Rat)
SwissProt	P49675(Human); P51557(Mouse); P97826(Rat)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****STAR Rabbit Polyclonal Antibody**

CAT. NO. APA13720

**DATA**

Western blot analysis of STAR expression in Jurkat (A), PC12 (B), mouse kidney (C) whole cell lysates. (Predicted band size: 31 kD; Observed band size: 31 kD)

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.