

DATASHEET

CD97 Rabbit Polyclonal Antibody

CAT. NO. APA14311

KEY FEATURES

Target	CD97	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

Receptor potentially involved in both adhesion and signaling processes early after leukocyte activation. Plays an essential role in leukocyte migration.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to CD97
Specificity	Recognizes endogenous levels of CD97 protein.
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human CD97
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 81; Observed: 105 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	CD97 antigen; Leukocyte antigen CD97; CD97
Gene Symbol	CD97
Entrez Gene	976(Human); 26364(Mouse)
SwissProt	P48960(Human); Q9Z0M6(Mouse)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**CD97 Rabbit Polyclonal Antibody**

CAT. NO. APA14311

DATA

Western blot analysis of CD97 expression in Jurkat (A), mouse lung (B), rat liver (C) whole cell lysates. (Predicted band size: 81; 86; 91 kD; Observed band size: 105 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.