

DATASHEET

Tropomodulin-3 Rabbit Polyclonal Antibody

CAT. NO. APA14407

KEY FEATURES

Target	Tropomodulin-3	Source / Host	Rabbit
Reactivity	Human, Mouse	Clonality	Polyclonal
Applications	WB, IHC	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.	Storage	at-20°C

BACKGROUND

Blocks the elongation and depolymerization of the actin filaments at the pointed end. The Tmod/TM complex contributes to the formation of the short actin protofilament, which in turn defines the geometry of the membrane skeleton .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:1000
IHC	1:50 - 1:200

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

Description	Rabbit polyclonal antibody to Tropomodulin-3
Specificity	Recognizes endogenous levels of Tropomodulin-3 protein.
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human Tropomodulin-3
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 39 kD; Observed: 40 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	Tropomodulin-3; Ubiquitous tropomodulin; U-Tmod
Gene Symbol	TMOD3
Entrez Gene	29766(Human); 50875(Mouse)
SwissProt	Q9NYL9(Human); Q9JHJ0(Mouse)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.


*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET

Tropomodulin-3 Rabbit Polyclonal Antibody

CAT. NO. APA14407

DATA



Western blot analysis of Tropomodulin-3 expression in A549 (A), HeLa (B), Jurkat (C) whole cell lysates. (Predicted band size: 39 kD; Observed band size: 40 kD)

Data 2

Immunohistochemical analysis of Tropomodulin-3 staining in human brain formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.