

DATASHEET

ARHGEF11 Rabbit Polyclonal Antibody

CAT. NO. APA14715

KEY FEATURES

Target	ARHGEF11	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB, IHC	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at-20°C

BACKGROUND

May play a role in the regulation of RhoA GTPase by guanine nucleotide-binding alpha-12 (GNA12) and alpha-13 (GNA13). Acts as guanine nucleotide exchange factor (GEF) for RhoA GTPase and may act as GTPase-activating protein (GAP) for GNA12 and GNA13. Involved in neurotrophin-induced neurite outgrowth.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
IHC	1:50 - 1:200

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

Description	Rabbit polyclonal antibody to ARHGEF11
Specificity	Recognizes endogenous levels of ARHGEF11 protein.
Antibody Type	Primary antibody
Immunogen	KLH-conjugated synthetic peptide of human ARHGEF11
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 167; Observed: 168 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	KIAA0380; Rho guanine nucleotide exchange factor 11; PDZ-RhoGEF
Gene Symbol	ARHGEF11
Entrez Gene	9826(Human); 78966(Rat)
SwissProt	O15085(Human); Q9ES67(Rat)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.


*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET

ARHGEF11 Rabbit Polyclonal Antibody

CAT. NO. APA14715

DATA



Western blot analysis of ARHGEF11 expression in mouse brain (A) whole cell lysates. (Predicted band size: 167; 172 kD; Observed band size: 168 kD)

Data 2

Immunohistochemical analysis of ARHGEF11 staining in human lung cancer formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.