

DATASHEET

USP39 Rabbit Polyclonal Antibody

CAT. NO. APA14891

KEY FEATURES

Target	USP39	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

Deubiquitinating enzyme that plays a role in many cellular processes including cellular antiviral response, epithelial morphogenesis, DNA repair or B-cell development . Plays a role in pre-mRNA splicing as a component of the U4/U6-U5 tri-snRNP, one of the building blocks of the precatalytic spliceosome . Specifically regulates immunoglobulin gene rearrangement in a spliceosome-dependent manner, which involves modulating chromatin interactions at the Igh locus and therefore plays an essential role in B-cell development . Regulates AURKB mRNA levels, and thereby plays a role in cytokinesis and in the spindle checkpoint . Regulates apoptosis and G2/M cell cycle checkpoint in response to DNA damage by deubiquitinating and stabilizing CHK2 .

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
----	----------------

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW


Description	Rabbit polyclonal antibody to USP39
Specificity	Recognizes endogenous levels of USP39 protein.
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human USP39
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 53; Observed: 65 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	U4/U6.U5 tri-snRNP-associated protein 2; Inactive ubiquitin-specific peptidase 39; SAD1 homolog; U4/U6.U5 tri-snRNP-associated 65 kDa protein; 65K
Gene Symbol	USP39
Entrez Gene	10713(Human); 28035(Mouse)
SwissProt	Q53GS9(Human); Q3TIX9(Mouse)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arex.bio or your local distributor.

*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET**USP39 Rabbit Polyclonal Antibody**

CAT. NO. APA14891

DATA

Western blot analysis of USP39 expression in HeLa (A), HT29 (B), mouse spleen (C), rat testis (D) whole cell lysates. (Predicted band size: 53; 56; 65 kD; Observed band size: 65 kD)

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.