

**DATASHEET**

**UACA Rabbit Polyclonal Antibody**

CAT. NO. APA14896

**KEY FEATURES**

Target	UACA	Source / Host	Rabbit
Reactivity	Human, Mouse	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

**BACKGROUND**

Regulates APAF1 expression and plays an important role in the regulation of stress-induced apoptosis. Promotes apoptosis by regulating three pathways, apoptosome up-regulation, LGALS3/galectin-3 down-regulation and NF-kappa-B inactivation. Regulates the redistribution of APAF1 into the nucleus after proapoptotic stress. Down-regulates the expression of LGALS3 by inhibiting NFKB1 .

**APPLICATION**

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
----	----------------

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

**PRODUCT OVERVIEW**


Description	Rabbit polyclonal antibody to UACA
Specificity	Recognizes endogenous levels of UACA protein.
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human UACA
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 161; Observed: 162 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	KIAA1561; Uveal autoantigen with coiled-coil domains and ankyrin repeats
Gene Symbol	UACA
Entrez Gene	55075(Human); 72565(Mouse)
SwissProt	Q9BZF9(Human); Q8CGB3(Mouse)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****UACA Rabbit Polyclonal Antibody**

CAT. NO. APA14896

**DATA**

Western blot analysis of UACA expression in HEK293T (A), SKOV3 (B) whole cell lysates.  
(Predicted band size: 161; 162 kD; Observed band size: 162 kD)

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.