

**DATASHEET**

**ADH1 Rabbit Polyclonal Antibody**

CAT. NO. APA15026

**KEY FEATURES**

Target	ADH1	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

**BACKGROUND**

Alcohol dehydrogenase . Oxidizes primary as well as secondary alcohols. Ethanol is a very poor substrate .

**APPLICATION**

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
----	----------------

\*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

**PRODUCT OVERVIEW**


Description	Rabbit polyclonal antibody to ADH1
Specificity	Recognizes endogenous levels of ADH1 protein
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human ADH1. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 39 kD; Observed: 40 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	ADH1; Alcohol dehydrogenase 1A; Alcohol dehydrogenase subunit alpha
Gene Symbol	ADH1A
Entrez Gene	124(Human)
SwissProt	P07327(Human)

\*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact [info@arexbio.com](mailto:info@arexbio.com) or your local distributor.

\*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

**DATASHEET****ADH1 Rabbit Polyclonal Antibody**

CAT. NO. APA15026

**DATA**

Western blot analysis of ADH1 expression in LO2 (A), rat liver (B) whole cell lysates.  
(Predicted band size: 39 kD; Observed band size: 40 kD)

**STORAGE**

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

**NOTE**

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.