

DATASHEET

TIP120B Rabbit Polyclonal Antibody

CAT. NO. APA15208

KEY FEATURES

Target	TIP120B	Source / Host	Rabbit
Reactivity	Human, Mouse, Rat	Clonality	Polyclonal
Applications	WB, IHC	Conjugation	Unconjugated
Form / Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.		Storage at -20°C

BACKGROUND

Probable assembly factor of SCF (SKP1-CUL1-F-box protein) E3 ubiquitin ligase complexes that promotes the exchange of the substrate-recognition F-box subunit in SCF complexes, thereby playing a key role in the cellular repertoire of SCF complexes.

APPLICATION

To ensure optimal assay performance, AREX recommends conducting reagent titration tailored to each testing system for optimal detection results.

WB	1:500 - 1:2000
IHC	1:50 - 1:200

*Results are sample-specific. Please refer to your local assay conditions and test parameters for reference.

PRODUCT OVERVIEW

Description	Rabbit polyclonal antibody to TIP120B
Specificity	Recognizes endogenous levels of TIP120B protein
Antibody Type	Primary antibody
Immunogen	Recombinant fusion protein of human TIP120B. The exact sequence is proprietary.
Purification	The antibody was purified by immunogen affinity chromatography.
Molecular Weight	Predicted: 122; Observed: 130 kD
Form/Buffer	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Alternative Names	KIAA0667; TIP120B; Cullin-associated NEDD8-dissociated protein 2; Cullin-associated and neddylation-dissociated protein 2; Epididymis tissue protein Li 169; TBP-interacting protein of 120 kDa B; TBP-interacting protein 120B; p120 CAND2
Gene Symbol	CAND2
Entrez Gene	23066(Human); 67088(Mouse); 192226(Rat)
SwissProt	O75155(Human); Q6ZQ73(Mouse); Q9R0L4(Rat)

*AREX continuously optimizes our products. Webpage content may not reflect the latest updates. For inquiries, please contact info@arexbio.com or your local distributor.


*Clone Number, Reactivity, Source/Host and Clonality can be found in the product name and Key Features section above.

DATASHEET


TIP120B Rabbit Polyclonal Antibody

CAT. NO. APA15208

DATA



Western blot analysis of TIP120B expression in Jurkat (A), NIH3T3 (B), mouse heart (C), rat heart (D) whole cell lysates. (Predicted band size: 122; 135 kD; Observed band size: 130 kD)



Immunohistochemical analysis of TIP120B staining in mouse spleen formalin fixed paraffin embedded tissue section. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH 6.0). The section was then incubated with the antibody at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

STORAGE

Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

NOTE

For Research Use Only. Not for diagnostic, therapeutics, prophylactic or in vivo use.